Locomotor Activity Test IACUC Standard Procedure Effective Date: October 2021



Description of Procedure:

The locomotor activity test (open field) is a simple means of assessing spontaneous locomotor activity and arousal in mice, and is necessary to establish baseline levels of motor activity.

Procedure Steps:

- 1. Mice are brought into the testing room under normal lighting conditions and allowed up to 30- 60 minutes of acclimation.
- 2. The chamber is cleaned with Vimoba prior to use.
- 3. The animals are placed into the locomotor activity chamber with infrared beams crossing the X, Y and Z axes that plot their ambulatory and fine motor movements, as well as rearing behavior. The animals are under no motivational constraints and are free to explore the environment. Alternatively, this test can be done with any open field enclosure, and the animal's movements and behavior can be video-recorded using a camera placed overhead. The researchers can then use post-video processing software for assessment of relevant variables.
- 4. The chamber is cleaned with 70% ethanol in between the testing of each mouse.