

Description of procedure:

In order to assess limb strength and coordination, a mouse is trained to hang by its forelimbs suspended from a wire or small bar.

Supplies:

- 1-2mm metal bar or wire
- Padding made from material that is easy to maintain/disinfect

Procedure Steps:

- 1. Padding is positioned under the bar so that the animal is not harmed upon releasing grip on the bar.
- 2. The mouse is conditioned to hang by its forelimbs suspended from a wire or small bar. The time that the animal maintains a grip on the wire or small bar is recorded, not to exceed 120 seconds.